

LANA Anaesthetic Machine

“Anaesthesia for all animals – LARGE and small”

Customers

- Vets
- Zoological parks
- Farms

Product Feature	Customer Benefit
Two machines in one	<ul style="list-style-type: none"> • Enables the administration of anesthesia to LARGE and small animals from the one machine • Convenient • Cost effective
Carbon dioxide absorber (180° rotatable)	<ul style="list-style-type: none"> • Ease of emptying and refilling • Need to lift bulky absorbers is eliminated
Hose-in-hose breathing circuit	<ul style="list-style-type: none"> • Less resistance than conventional breathing circuits • Lower dead space and pressure tested to prevent leaks • Significantly lower operating costs
KDK™ flow meter	<ul style="list-style-type: none"> • Variable oxygen flow rates – between ½ to 8 L/min • Oxygen flush button to quickly fill face mask reservoir bag with oxygen • Superior quality - all high pressure parts made from brass
Breathing monitor	<ul style="list-style-type: none"> • Monitors the respiration of the patient, and produces an audible beep with each breath • Security and peace of mind
Trolley mounted	<ul style="list-style-type: none"> • Can be easily transported and used in any location • Includes support for medical oxygen cylinder

Key Innovations

- Combined anesthetic machine for LARGE and small animals
- Cost effective (i.e. two machines in one)
- Ease of emptying and refilling carbon dioxide absorber (under two minutes)

Summary

- An innovative medical device that enables the administration of anesthesia to LARGE and small animals from the one machine