

KAB™ Carbon Dioxide Absorber

“Enhance your breathing circuit”



Recommended

- For the provision of a clean and fresh breathing circuit for every patient

Key Product Highlights

- Easily removed from breathing circuit
 - Enables rapid changeover and replacement (under one minute)
 - Prevents cross infection between patients
 - Disposable version / Re-usable version
 - Single patient use / Simple sterilisation by autoclave
 - Pre-filled with carbon dioxide absorbing agent / Refillable via plug at back of device
 - Low resistance inspiratory and expiratory valves – made from medical grade silicone
 - Automatic exhaust valve – APL valve
 - Superior quality – made from medical grade polycarbonate
 - Features international standard fittings
- CE marked, FDA # K953411, US patent # 5,556,669 and Australia patent # 694,208

Usage and Directions

- Ensure that KAB™ CO₂ absorber is filled with CO₂ absorbing agent (e.g. soda lime)
- Attach to breathing circuit and sterilise or dispose of between patients