

Tourniquet

“Quick and easy”

Customers

- Hospitals
- Pathology
- General Practitioners (e.g. Doctors)

Product Feature	Customer Benefit
Seat belt design	<ul style="list-style-type: none"> • Click in and out • Enables quick application and release
Latex free	<ul style="list-style-type: none"> • Non-allergenic
Easy fitting	<ul style="list-style-type: none"> • Relatively comfortable for the patient • Easy to apply by the operator
Rapid release lever	<ul style="list-style-type: none"> • Enables for quick release
Slow release lever	<ul style="list-style-type: none"> • Enables gradual and measured release
Made from medical grade polycarbonate	<ul style="list-style-type: none"> • Robust and long-lasting • Can be sterilised by autoclaving

Key Innovations

- Latex free and non-allergenic
- Seat belt design with slow and rapid release levers
- Can be easily sterilised by autoclaving

Summary

- An ingenious and simple-to-use medical device for the taking of blood from a patient