

# Tourniquet

“Quick and easy”



## Recommended

- For the obstruction of venous flow to enable the taking of blood from a patient

## Key Product Highlights

- Easy fitting tourniquet
- Seat belt design for instant application and release
- Latex free and non-allergenic
- Slow and rapid release levers
- Robust – made from medical grade polycarbonate
- Can be easily sterilised by autoclaving
- Suitable for children and adults of all ages
- Can be used as a gentle restraint for limbs during transport, resuscitation or anaesthesia

## Usage and Directions

- Apply tourniquet to area where blood is required
- Adjust tension using levers, and release when blood has been taken