

BioMate®

“For the cleanest punch, everytime”



Recommended

- For the precise removal of clean biopsies and specimens

Key Product Highlights

- Razor sharp stainless steel cutting blade
- Ribbed handle for maximum control and comfort
- Rigid plastic packaging for greater protection
- Hollow handle for the easy removal of biopsies and specimens
- Single use / disposable
- Wide range of sizes available
 - Sizes 2, 3, 4, 5, 6 and 8 mm

Usage and Directions

- Clean and / or anaesthetise area for biopsy removal
- Hold tissue firmly, stretch skin and push BioMate™ downward while rotating
- Remove specimen from BioMate™ for diagnosis or disposal
- Close skin opening in patient